



**Buy 2 cases of DURAN® YOUTILITY Bottles
Receive a FREE pack of Colored Bottle Tags**



DURAN® YOUTILITY Laboratory Bottle

The specially shaped gripping zones on both sides of the bottle enable easier and safer handling. With the new YOUTILITY bottle thread opening or closing the bottle is significantly faster. The thread is fully compatible with DIN GL 45 closures and other accessories. The slimmer YOUTILITY bottle shape allows a more optimal use of limited space in autoclaves and laboratory refrigerators. A pre-defined labelling area is compatible with the dedicated DURAN® self-adhesive YOUTILITY bottle labels. Nominal volume is shown at the top of the easy-to-read graduation scale for fast determination of the volumes. The glass bottle body is moulded from the tried and tested DURAN® borosilicate 3.3 pharmacopoeial Type 1 neutral glass. DURAN® glass offers very good chemical resistance and high temperature resistance. Each DURAN® YOUTILITY bottle is supplied as a complete system, with a pouring ring (PP¹) and a GL 45 cap (PP¹).

Cat. No. Clear Glass	Cat. No. Amber Glass	Capacity (mL)	Thread (GL)	d (mm)	h (mm)	Qty/ Case
12-21 881 2854	6-21 886 2859	125	45	55	124	4
12-21 881 3653	6-21 886 3658	250	45	66	158	4
12-21 881 4452	6-21 886 4457	500	45	78	193	4
12-21 881 5457	6-21 886 5453	1000	45	93	253	4

DURAN® Bottle Tag

The colourful silicone GL 45 Bottle Tags can be easily and securely attached around the neck of the YOUTILITY bottle. The Bottle Tags are available in several distinct colours for easy customization and unambiguous bottle identification. The elastic Bottle Tags will fit around any DIN GL 45 bottle neck, and can even be applied with the screw cap in place.

Cat. No.	Color	Qty/Case
12-29 243 2904	Eight colours. Two silicone Bottle Tags are supplied in each colour	16

To claim your redemption (redeem by 7/31/19)

1. Obtain your Proof of Purchase (dated between 1/1/19 - 6/30/19). Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit www.dwkdeals.com
3. Attach your Proof of Purchase and fill-out the online form.
4. Enter the part numbers for freeware. Free items must be of equal or lesser value. Visit www.DWK.com for more details. Not for resale. Limited to two redemptions per promotion, per customer and/or per laboratory.