



Save 15% on Tacta® and Proline® Pipette Multipacks

Top Features of the Tacta® Mechanical Pipette

- Superior Ergonomics
- Single-Handed Volume Locking Mechanism
- Large, Visible Display
- Easy Calibration for a Variety of Liquids

Popular Applications for the Tacta® Mechanical Pipette

- Pharmaceutical QA/QC
- Biopharma/Biotech R&D
- Academic Research

Top Features of Proline® Plus Mechanical Pipettes

- Comfortable, reliable liquid handling
- Easy to maintain and clean
- Fully autoclavable

Popular Applications for Proline® Plus Mechanical Pipettes

- Academic
- Industrial
- Municipal

Order Information

Sartorius Part No.	Description	List Price	Special Price
109-LH-729670	Tacta 3-pack 10	\$1,233.12	\$1,048.15
109-LH-729671	Tacta 3-pack 20	\$1,233.12	\$1,048.15
109-LH-729672	Tacta 4-pack 20	\$1,530.39	\$1,300.83
109-LH-729673	Tacta 4-pack 100	\$1,530.39	\$1,300.83
109-LH-729674	Tacta 5-pack 10	\$1,794.63	\$1,525.44
109-LH-729675	Tacta 5-pack 20	\$1,794.63	\$1,525.44
109-LH-729676	Tacta 3+1-pack	\$2,080.89	\$1,768.76
109-LH-729677	Tacta 3-pack Max	\$1,233.12	\$1,048.15

Sartorius Part No.	Description	List Price	Special Price
109-LH-728670	Proline Plus 3-pack 10	\$845.57	\$718.73
109-LH-728671	Proline Plus 3-pack 20	\$845.57	\$718.73
109-LH-728672	Proline Plus 4-pack 20	\$1,065.77	\$905.90
109-LH-728673	Proline Plus 4-pack 100	\$1,065.77	\$905.90
109-LH-728674	Proline Plus 5-pack 10	\$1,376.25	\$1,169.81
109-LH-728675	Proline Plus 5-pack 20	\$1,376.25	\$1,169.81
109-LH-728676	Proline Plus 3+1-pack	\$1,519.38	\$1,291.47
109-LH-728677	Proline Plus 3-pack Max	\$845.57	\$718.73

Feel the Comfort

Tacta® rests lightly in your hand thanks to its ergonomically designed handle and finger hook. The pipette is effortless to use and exceptionally comfortable to hold. They eases your workload and protects you from strain, even when you pipette for extended periods of your working day.

Volume Adjustment


Use the Sartorius Optilock® feature to prevent accidental volume changes during pipetting. Either hold down the volume lock button with one hand and adjust the volume with the other, then release to lock, or with just one hand, slide the volume lock up, adjust the volume and slide down to lock again.

Large, Clear Display

The large, easy-to-read display helps you to see all four digits of the set volume from a distance without straining your eyes. The volume is easy to read even when the pipette is angled – eliminating the need to turn your head into an uncomfortable position.

Calibration Adjustment


Tacta® is very easy to adjust, for various types of liquids, using a simple adjustment key. The adjustment scale shows the degree of adjustment. By noting this value for a specific liquid, you can return to that setting any time. Easy calibration adjustment ensures accurate results with a variety of liquids.



Optilock®



The volume is easy to read even when the pipette is angled.



Calibration Adjustment

Special offer valid only for the US until December 31, 2026 and cannot be combined with other offers.

Contact:

www.krackeler.com
800-334-7725
sales@krackeler.com

